



|     |      |       |       |             |  |   |  |  |
|-----|------|-------|-------|-------------|--|---|--|--|
| 1   |      |       | 2     |             |  | 3 |  |  |
| MOD | DATE | DRND. | APPD. | DESCRIPTION |  |   |  |  |





REMOVE ALL BURRS AND SHARP EDGES.

|                       |  |          |           |  |         |
|-----------------------|--|----------|-----------|--|---------|
|                       |  | PEEK     |           |  | 0,00018 |
| DESCRIPTION / REMARKS |  | MATERIAL | TREATMENT |  | W.(Kg)  |

SCALE: 5/ 1


0mm

20

|       |        |          |
|-------|--------|----------|
|       | NAME   | DATE     |
| DRN.  | HEYMAN | 27/07/17 |
| CKD.  |        |          |
| APPD. |        |          |

BM23 - EXAFS DETECTORS  
CRYSTAL ANALYSER-CRYOSTAT  
SAMPLE HOLDER

GEN. TOL. : ISO2768 f H  
GEN. SURFACE FINISH : Ra 1,6



EUROPEAN SYNCHROTRON  
RADIATION FACILITY

CS 40220 - 38043 Grenoble Cedex 9  
TEL. 04-76-88-20-00 WWW.ESRF.FR

A4

ASSY. 61.56.1074

61.56.1101

|   |   |              |
|---|---|--------------|
| 1 | 2 | 3 SolidWorks |
|---|---|--------------|